

SANI SYSTEM POLTI "CHECK"



12 HP MED 50ml single-use bottles
(included in the packaging)

Main features

- Electronic system for all functions' management and control;
- 4 preset sanitization programs;
- User identification system through RFID cards;
- Built-in thermal printer;
- Numeric interface with back-lit LCD display;
- Steam nozzle with worldwide Polti patented controlled expansion chamber to achieve a high temperature steam outflow;
- 180° overheated dry steam outflow;
- Max steam release: 100 g/min;
- HP MED minimum guaranteed consumption: 0,4 ml/min;
- Max pressure: 6 bar;
- Pressure gauge;
- Boiler with auto-filling system for unlimited operating time;
- Reinforced Aisi 316 stainless steel boiler;
- Connection for HP MED sanitizing solution single-use bottle placed under the steam nozzle;
- Manual cable winder on the back of the machine;
- Stainless steel handle for easy transport;
- 4 anti-trace revolving wheels (1 antistatic) with brakes;

- 12 HpMed single-use bottles (50ml each) included in the packaging;
- Adapter for a concentrated jet of steam;



Spacer



Adapter for a concentrated steam jet

Technical Data

- Voltage: 220-240 V ~ 50-60 Hz
- Flex hose length: 2,5 m
- Water tank capacity: 5 l
- Boiler power: 2200 W
- Steam nozzle power: 1500 W
- Total power: Max 2250 W
- Net weight: 29 Kg
- Dimensions (LxWxH): 47 x 45 x 107 cm
- Electric cable length: 5 m
- Certification: Class IIA CE Medical

	COD. PRODOTTO	COD. EAN 13
SANI SYSTEM CHECK	PTEU0207	8007411006784
HPMED (12x50ml)	PAEU0177	8007411802324