

# SANI SYSTEM POLTI



12 HpMed 50ml single-use bottles  
(included in the packaging)

## Main Features

- Steam nozzle with worldwide Polti patented controlled expansion chamber to achieve a high temperature steam outflow;
- 180° overheated dry steam outflow;
- Continuous steam quantity adjustment from 0 to 100 g/min;
- HpMed minimum guaranteed consumption: 0,4 ml/min;
- Max pressure: 6 bar;
- Pressure gauge;
- Boiler with auto-filling system for unlimited operating time;
- Reinforced Aisi 316 stainless steel boiler;
- Connection for HPMED sanitizing solution single-use bottle placed under the steam nozzle;
- Manual cable winder on the back of the machine;
- Stainless steel handle for easy transport;
- 4 anti-trace revolving wheels (1 antistatic) with brakes;
- 12 HpMed single-use bottles (50ml each) included in the packaging;
- Adapter for a concentrated steam jet;



Spacer



Adapter for a  
concentrated steam jet

### **Technical Data**

- Voltage: 220-240 V ~ 50-60 Hz
- Flex hose length: 2,5 m
- Water tank capacity: 5 l
- Boiler power: 2200 W
- Steam nozzle power: 1500 W
- Total power: Max 2500 W
- Net weight: 27,5 kg
- Dimensions (LxWxH): 47 x 45 x 90,5 cm
- Electric cable length: 5 m
- Certification: Class IIA CE Medical

	<b>PRODUCT CODE</b>	<b>EAN 13 CODE</b>
<b>SANI SYSTEM 02</b>	<b>PTEU0205</b>	<b>8007411006586</b>
<b>HP MED (12x50ml)</b>	<b>PAEU0177</b>	<b>8007411802324</b>