



SANI SYSTEM EXPRESS



*12HPMED 50 ml disposable bottles
(included in the packaging)*

Main features

- Steam nozzle with worldwide Polti patented controller expansion chamber to achieve a high temperature steam out flow
- Steam release button with safety lock device .
- up to 180°C overheated dry steam out flow
- Max steam release : 100 g/min
- max. pressure 5 bar;
- Boiler with auto-filling system for unlimited operating time
- aluminium Extra Alp boiler
- Connection for HPMed sanitising solution bottle placed under the steam nozzle
- pull-out handle for easy transport
- Hooks for parking the tubes and the cable
- 4 anti-trace revolving wheels (1 antistatic)
- 12 HpMed disposable bottles (50 ml each) included in the packaging
- adapter for a concentrated jet of steam
- side pocket for storing accessories



MEDICAL DIVISION



Spacer



Adapter for a concentrated steam jet

Technical data

- Voltage: 220-240 V ~ 50-60 Hz
- Flex hose length 2,5 m
- Water tank capacity: 3 l
- Boiler power: 1000 W
- Steam nozzle power : 1500 W
- Total power: Max 2650 W
- Dimensions (DxH): 45 x 39 cm
- Height of transport handle (min / max): 49 cm / 68 cm
- Electric cable length 4 m
- Certification: Class IIA CE Medicale

	PRODUCT CODE	EAN COD.
SANI SYSTEM express	PTEU0221	8007411008009
HP MED (12x50ml)	PAEU0177	8007411802324

